

Colored Clips

Seamless Syringe Identification

Custom Engraved Clips: Secure and Durable

Accurately identify your syringes with our customizable engraved clips. Whether you need quick and easy volume identification, fluid labeling such as CH₃OH, or other essential information, our advanced engraving technology ensures clear and permanent labeling.

Simple Clips: Versatile and Reliable

Designed to simplify your syringe identification process. With an array of three colors to choose from, these clips allow you to quickly and easily distinguish between different syringe sizes, applications, or fluids, serving as visual cues for efficient identification.

Explore the Numerous Benefits

Durability: Crafted from high-quality PEEK material, known for its exceptional strength.

Chemical Resistance: Ideal for use with most laboratory solvents used in HPLC, except for concentrated nitric acid, concentrated sulfuric acid, methylene chloride, DMSO, and THF.

Reusable: Autoclavable and resistant to alkalis and detergents, the clips are easy to clean and reuse.

Permanent Markings: The engravings are permanent and resistant to fading, ensuring the information remains visible and intact even with repetitive use and cleaning.

Time-Saving and Error-Reducing: Quickly locate the desired syringe, minimizing the risk of mix-ups and errors.



Precision Engraved Clips: Reusable and Professional

CLIP CUSTOMIZATION OPTIONS

- Volume Measurements:** Engrave the clips with volume indicators such as "1 mL," "5 mL," or "10 mL" to easily identify syringes of different capacities.
- Substance Names:** Customize the clips with substance names like "CH₃OH" (methanol), "NaCl" (sodium chloride), or "H₂O" (water) to distinguish syringes containing specific substances.
- Safety Labels:** Engrave warning symbols or labels like "Caution" or "Biohazard" to alert users about potential risks associated with specific syringes.
- Identification Codes:** Assign unique identification codes or numbers to each clip to facilitate inventory management and traceability of syringes.

RECOMMENDED ENGRAVING CONFIGURATIONS

- 5 µL/ 10 µL/ 25 µL/ 100 µL/ 1 mL / 2.5 mL: Arial 4 – BOLD = 2 rows x 7characters

25 µL/ 50 µL/ 100 µL/ 250 µL / 500 µL/ 1 mL/ 5 mL / 10 mL: Arial 4 - BOLD = 2 rows x 7characters
- 5 µL/ 10 µL/ 25 µL/ 100 µL: Arial 6 – BOLD =1 row x 5 characters

5 mL/ 10 mL: Arial 6 – BOLD = 1 row x 10 characters

ORDERING CONFIGURATIONS

Size (Syringe Barrel Outer Diameter)	Simple Clip PN, 5/pk	Custom Engraved Clip PN, 5/pk	Compatible Syringe Series	Color
Ø6.6 mm	10164247	10164252	775, 701, 1701, 1702* slim, 1710* slim	<div><div></div><div></div><div></div></div>
Ø7.7 mm	10164248	10164253	702, 705, 710, 725, 750, 1702*, 1705, 1710*, 1725, 1750, 1001	<div><div></div><div></div><div></div></div>
Ø9.65 mm	10164249	10164254	1002 , 1001	<div><div></div><div></div><div></div></div>
Ø13.4 mm	10164250	10164255	1005	<div><div></div><div></div><div></div></div>
Ø17.4 mm	10164251	10164256	1010	<div><div></div><div></div><div></div></div>

*1702 and 1710 syringe barrel diameters vary. To ensure compatibility, confirm the outer diameter of your syringe barrel.

© 2023 Hamilton Company. All rights reserved. All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries. Lit. No. 10166558/00 – 08/2023



To find a representative in your area, please visit:
www.hamiltoncompany.com/contact
Email: contact.hce.ro@hamilton-ce.com

Hamilton Europe, Middle East, Asia & Africa
Hamilton Central Europe S.R.L.
Str. Hamilton nr. 2-4
307210 Giarmata, Romania
Tel: +40-356-635-055

