



MSPE

Unique
technology

SPE
EXTRA



Speed up
your sample preparation

- Simple and fast sample preparation
- 2 in 1 – micro-solid phase extraction and filtration in one step
- SPE manifold replaced by use of centrifuge filters

Shorter sample
preparation time



FAST

Simple
technique



EASY



CLEAN

Lower solvent
consumption

SPE and filtration
at once using 0.22 μ m
Nylon membrane




COMPACT



CLEAR

Centrifuging
multiple samples at once

General sample preparation workflow



New products

Item description	Note	p/n
Micro Spin SpeExtra MSPE HLB column, 0.22 µm Nylon membrane, 15 mg, 0.7 ml, 50 pcs		MSP-5877-AA07
Micro Spin SpeExtra MSPE C18 column, 0.22 µm Nylon membrane, 30 mg, 0.7 ml, 50 pcs		MSP-5876-AB07
Micro Spin SpeExtra MSPE C18-P column, 0.22 µm Nylon membrane, 30 mg, 0.7 ml, 50 pcs	Coming soon	MSP-5879-AB07
Micro Spin SpeExtra MSPE HLB column, 0.22 µm Nylon membrane, 30 mg, 0.7 ml, 50 pcs		MSP-5877-AB07
Micro Spin SpeExtra MSPE C18 column, 0.22 µm Nylon membrane, 50 mg, 0.7 ml, 50 pcs		MSP-5876-AC07
Micro Spin SpeExtra MSPE C18-P column, 0.22 µm Nylon membrane, 50 mg, 0.7 ml, 50 pcs	Coming soon	MSP-5879-AB07

